



INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification ⁷ : G01N 33/00	A1	(11) International Publication Number: WO 00/58724 (43) International Publication Date: 5 October 2000 (05.10.00)
---	----	--

(21) International Application Number: PCT/EP00/02887

(22) International Filing Date: 28 March 2000 (28.03.00)

(30) Priority Data:
9907249.8 29 March 1999 (29.03.99) GB(71) Applicant (for all designated States except US): ECOLE
POLYTECHNIQUE FEDERALE DE LAUSANNE
[CH/CH]; Laboratoire d'Electrochimie, Département de
Chimie-ICP III, CH-1015 Lausanne (CH).

(72) Inventors; and

(75) Inventors/Applicants (for US only): ROSSIER, Joël, S.
[CH/CH]; Ch. du Chamossaire 2, CH-1860 Aigle (CH).
REYMOND, Frédéric [CH/CH]; Ch. des Marionnettes 15,
CH-1093 La Conversion (CH). GIRAULT, Hubert, H.
[FR/CH]; CH-1088 Ropraz (CH).(74) Agent: MANATON, Ross, Timothy; J.Y. & G.W. John-
son, Kingsbourne House, 229-231 High Holborn, London
WC1V 7DP (GB).(81) Designated States: AE, AG, AL, AM, AT, AU, AZ, BA, BB,
BG, BR, BY, CA, CH, CN, CR, CU, CZ, DE, DK, DM,
DZ, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL,
IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU,
LV, MA, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT,
RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, TZ,
UA, UG, US, UZ, VN, YU, ZA, ZW, ARIPO patent (GH,
GM, KE, LS, MW, SD, SL, SZ, TZ, UG, ZW), Eurasian
patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European
patent (AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR,
IE, IT, LU, MC, NL, PT, SE), OAPI patent (BF, BJ, CF,
CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG).

Published


With international search report.

Before the expiration of the time limit for amending the
claims and to be republished in the event of the receipt of
amendments.


(54) Title: MICROSCALE TOTAL ANALYSIS SYSTEM

(57) Abstract

The invention relates to an apparatus for performing chemical assays in an aqueous medium. The apparatus contains a reaction chamber(s) and a liquid in-flow channel connected to each chamber. The flow of liquid through the fluid in-flow channel to the reaction chamber is controlled by the presence of a hydrophobic inner surface on the walls of the in-flow channel. Under normal conditions fluid will not flow through the channel. However, application of an external force pushes the liquid through said channel into the reaction chamber. The invention is applicable to the monitoring of many different molecular interactions, in particular molecular recognition between an immobilised affinity partner and a species in solution, such as immunoglobulin/antigen interaction, DNA hybridisation, haptamer-protein interaction, drug and virus detection, high throughput screening of synthetic molecules and for determining the concentration and reaction kinetics of target species.



A



B